



Recombinant Enterobacteria phage T4 Single-stranded DNA-binding protein (32)

Product Code	CSB-BP022706EDZ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P03695
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage T4 (Bacteriophage T4)
Purity	≥85% (SDS-PAGE)
Sequence	MFKRKSTAE L AAQMAKLNGN KGFSSSEDKGE WKLKLDNAGN GQAVIRFLPS KNDEQAPFAI LVNHGFKKNG KWIYETCSST HGDYDSCPVC QYISKNDLYN TDNKEYSLVK RKTSYWANIL VVKDPAAPEN EGKVFYRFG KKIWDKINAM IAVDVEMGET PVDVTCPWEG ANFVLKVKQV SGFSNYDESK FLNQSAIPNI DDES FQKELF EQMVDLSEMT SKDKFKSFEE LNTKFGQVMG TAVMGGAAAT AAKKADKVAD DLDAFNVDDF NTKTEDDFMS SSSGSSSSAD DTDLDDLLND L
Source	Baculovirus
Target Names	32
Protein Names	Recommended name: Single-stranded DNA-binding protein Alternative name(s): Gp32 Helix-destabilizing protein
Expression Region	1-301
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.