



# Recombinant Human Alpha-1B-glycoprotein (A1BG)

<b>Product Code</b>	CSB-EP001001HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P04217
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AIFYETQPS LWAEESELLK PLANVTLTCQ AHLETPDFQL FKNGVAQEPV HLDSPAIKHQ FLLTGDTQGR YRCRSLSTG WTQLSKLLEL TGPKSLPAPW LSMAPVSWIT PGLKTTAVCR GVLRGVTFL RREGDHEFLE VPEAQEDVEA TFPVHQPGNY SCSYRTDGEG ALSEPSATVT IEELAAPPPP VLMHHGESSQ VLHPGNKVTL TCVAPLSGVD FQLRRGEKEL LVPRSSTSPD RIFFHLNAVA LGDGGHYTCR YRLHDNQNGW SGDSAPVELI LSEDTLPAPE FSPEPESGRA LRLRCLAPLE GARFALVRED RGGRRVHRFQ SPAGTEALFE LHNISVADSA NYSCVYVDLK PPFGGSSAPSE RLELHVDGPP PRPQLRATWS GAVLAGRDAV LRCEGPIPDV TFELLREGET KAVKTVRTPG AAANLELIFV GPQHAGNYRC RYRSWVPHTF ESELSDPVEL LVAES
<b>Source</b>	E.coli
<b>Target Names</b>	A1BG
<b>Protein Names</b>	Recommended name: Alpha-1B-glycoprotein Alternative name(s): Alpha-1-B glycoprotein
<b>Expression Region</b>	22-495
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.