



Recombinant Human Alpha-1B-glycoprotein (A1BG)

Product Code	CSB-YP001001HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P04217
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AIFYETQPS LWAESSELLK PLANVTLTCC AHLETPDFQL FKNGVAQEPV HLDSIPAICKHQ FLLTGDTQGR YRCRSLSTG WTQLSKLLEL TGPKSLPAPW LSMAPVSWIT PGLKTTAVCR GVLRGVTFL RREGDHEFLE VPEAQEDVEA TFPVHQPGNY SCSYRTDGEG ALSEPSATVT IEELAAPPPP VLMHHGESSQ VLHPGNKVTL TCVAPLSGVD FQLRRGEKEL LVPRSSTSPD RIFFHLNAVA LGDGGHYTCR YRLHDNQNGW SGDSAPVELI LSETLPAPPE FSPEPESGRA LRLRCLAPLE GARFALVRED RGGRRVHRFQ SPAGTEALFE LHNISVADSA NYSCVYVDLK PPFGGSSAPSE RLELHVDGPP PRPQLRATWS GAVLAGRDAV LRCEGPIPDV TFELLREGET KAVKTVRTPG AAANLELIFV GPQHAGNYRC RYRSWVPHTF ESELSDPVEL LVAES
Source	Yeast
Target Names	A1BG
Protein Names	Recommended name: Alpha-1B-glycoprotein Alternative name(s): Alpha-1-B glycoprotein
Expression Region	22-495
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.