



Recombinant Human Protamine-2 (PRM2)

Product Code	CSB-EP018729HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P04554
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	VRYSRVSLS ERSHEVYRQQ LHGQEQGHHG QEEQGLSPEH VEVYERTHGQ SHYRRRHCSR RRLHRIHRRQ HRSCRRRKRR SCRHHRRHR GCRTRKRTCR RH
Source	E.coli
Target Names	PRM2
Protein Names	Recommended name: Protamine-2 Alternative name(s): Sperm histone P2 Sperm protamine P2 Cleaved into the following 7 chains: 1. Basic nuclear protein HPI1 2. Basic nuclear protein HPI2 3. Basic nuclear protein HPS1 4. Basic nucl
Expression Region	2-102
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.