



# Recombinant Human Plasma serine protease inhibitor (SERPINA5)

<b>Product Code</b>	CSB-BP021061HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P05154
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	EMKKR VEDLHVGATV APSSRRDFTF DLYRALASAA PSQSIFFPV SISMSLAMLS LGAGSSTKMQ ILEGLGLNLQ KSSEKELHRG FQQLQELNQ PRDGFQLSLG NALFTDLVVD LQDTFVSAMK TLYLADTFPT NFRDSAGAMK QINDYVAKQT KGKIVDLLKN LDSNAVIMV NYIFFKAKWE TSFNHKGTE QDFYVTSETV VRVPMMSRED QYHYLLDRNL SCRVVGVPYQ GNATALFILP SEGKMQQVEN GLSEKTLRKW LKMFKKRQLE LYLPKFSIEG SYQLEKVLPS LGISNVFTSH ADLSGISNHS NIQVSEMVHK AVVEVDESGT RAAATGTIF TFRSARLNSQ RLVFNRPFLM FIVDNNILFL GKVNRP
<b>Source</b>	Baculovirus
<b>Target Names</b>	SERPINA5
<b>Protein Names</b>	Recommended name: Plasma serine protease inhibitor Alternative name(s): Acrosomal serine protease inhibitor Plasminogen activator inhibitor 3 Short name= PAI-3 Short name= PAI3 Protein C inhibitor Short name= PCI S
<b>Expression Region</b>	26-406
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.