



Recombinant Mouse Transcription factor SOX-2 (Sox2)

Product Code	CSB-EP022426MO-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P48432
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MYNMMETELK PPGPQQASGG GGGGGNATAA ATGGNQKNSP DRVKRPMNAF MVWSRGQRRK MAQENPKMHN SEISKRLGAE WKLLSETEKR PFIDEAKRLR ALHMKEHPDY KYRPRRKT KT LMKKDKYTLP GGLLAPGGNS MASGVGVGAG LGGGLNQRMDSYAHMNGWSN GSYSMMQEQL GYPQHPLNA HGAAQMMPMH RYVVSALQYN SMTSSQTYMN GSPTYMSYS QQGTPGMALG SMGSVVKSEA SSSPPVVTSS SHSRAPCQAG DLRDMISMYL PGAEVPEPAA PSRLHMAQHY QSGPVPGTAK YGTLPLSHM</p>
Source	E.coli
Target Names	Sox2
Protein Names	Recommended name: Transcription factor SOX-2
Expression Region	1-319
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This intronless gene encodes a member of the SRY-related HMG-box (SOX) family of transcription factors involved in the regulation of embryonic development and in the determination of cell fate. The product of this gene is required for stem-cell maintenance in the central nervous system, and also regulates gene expression in the stomach. Mutations in this gene have been associated with optic nerve hypoplasia and with syndromic microphthalmia, a severe form of structural eye malformation. This gene lies within an intron of another gene called SOX2 overlapping transcript (SOX2OT).
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.