



Recombinant Acyl carrier protein (acpP)

Product Code	CSB-EP015636EGX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0A6A9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	STIEERVKK IIGEQLGVKQ EEVTNNASFV EDLGADSLDT VELVMALEEE FDTEIPDEEA EKITTVQAAI DYINGHQA
Source	E.coli
Target Names	acpP
Protein Names	Recommended name: Acyl carrier protein Short name= ACP Alternative name(s): Cytosolic-activating factor Short name= CAF Fatty acid synthase acyl carrier protein
Expression Region	2-78
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.