




Recombinant Escherichia coli O157:H7 Peptide deformylase (def)

Product Code	CSB-EP017707EOD
Relevance	Removes the formyl group from the N-terminal Met of newly synthesized proteins. Requires at least a dipeptide for an efficient rate of reaction. N-terminal L-methionine is a prerequisite for activity but the enzyme has broad specificity at other positions.
Abbreviation	Recombinant E.coli O157:H7 def protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0A6K5
Alias	Polypeptide deformylase
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	SVLQVLHIPDERLRKVKAKPVVEEVNAEIQRIVDDMFETMYAEEGIGLAATQVDIH QRIIVIDVSENRLVRLINPELLEKSGETGIEEGCLSIPEQRALVPRAEKVKIRA LDRDGKPFLEADGLLAICIQHMDHLVGLFMDYLSPLKQQRIRQKVEKLDRL KARA
Source	E.coli
Target Names	def
Expression Region	2-169aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	35.2kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.