



# Recombinant Human Glutathione synthetase (GSS)

<b>Product Code</b>	CSB-EP009969HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P48637
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	<p>ATNWGSLQ DKQQLEELAR QAVDRALAEG VLLRTSQEPT SSEVVSYPF            TLFPSLVPSA LLEQAYAVQM DFNLLVDAVS QNAAFLEQTL SSTIKQDDFT            ARLFDIHKQV LKEGIAQTVF LGLNRSYMF QRSADGSPAL KQIEINTISA            SFGGLASRTP AVHRHVLSVL SKTKEAGKIL SNNPSKGLAL GIAKAWELYG            SPNALVLLIA QEKERNIFDQ RAIENELLAR NIHVIRRTFE DISEKGLDQ            DRRLFVDGQE IAVVYFRDGY MPRQYSLQNW EARLLERSH AAKCPDIATQ            LAGTKKVQQE LSRPGMLEML LPGQPEAVAR LRATFAGLYS LDVGEEGDQA            IAEALAAPSR FVLKPQREGG GNNLYGEEMV QALKQLKDSE ERASYILMEK            IEPEPFENCL LRPGSPARVV QCISELGIFG VYVRQEKTLV MNKHVGHLLR            TKAIEHADGG VAAGVAVLND PYPV</p>
<b>Source</b>	E.coli
<b>Target Names</b>	GSS
<b>Protein Names</b>	Recommended name: Glutathione synthetase Short name= GSH synthetase Short name= GSH-S EC= 6.3.2.3 Alternative name(s): Glutathione synthase
<b>Expression Region</b>	2-474
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	Glutathione is important for a variety of biological functions, including protection of cells from oxidative damage by free radicals, detoxification of xenobiotics, and membrane transport. This protein functions as a homodimer to catalyze the second step of glutathione biosynthesis, which is the ATP-dependent conversion of gamma-L-glutamyl-L-cysteine to glutathione. Defects in this gene are a cause of glutathione synthetase deficiency.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.