



Recombinant Escherichia coli O6:H1 ATP synthase subunit alpha (atpA)

Product Code	CSB-BP002344EGX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0ABB1
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MQLNSTEISE LIKQRIAQFN VVSEAHNEGT IVSVSDGVIR IHGLADCMQG EMISLPGNRY AIALNLERDS VGAVVMGPYA DLAEGMKVKC TGRILEVPVG RGLLGRVVNT LGAPIDGKGP LDHDFSAVE AIAPGVIERQ SVDQPVTGY KAVDSMIPIG RGQRELIIGD RQTGKTALAI DAIINQRDSG IKCIYVAIGQ KASTISNVVR KLEEHEGALAN TIVVATASE SAALQYLAPY AGCAMGEYFR DRGEDALIY DDLISKQAVAY RQISLLLRP PGREAFPGDV FYLHSRLLER AARVNAEYVE AFTKGEVKGK TGSLTALPII ETQAGDVSAF VPTNVISITD GQIFLETNLF NAGIRPAVNP GISVSRVGGGA AQTKIMKKLS GGIRTALAQY RELAAFSQFA SDLDDATRKQ LDHGQKVTEL LKQKQYAPMS VAQQSLVLFA AERGLADVE LSKIGSFEAA LLAYVDRDHA PLMQEINQTG GYNDEIEGKL KGILDSFKAT QSW
Source	Baculovirus
Target Names	atpA
Protein Names	Recommended name: ATP synthase subunit alpha EC= 3.6.3.14 Alternative name(s): ATP synthase F1 sector subunit alpha F-ATPase subunit alpha
Expression Region	1-513
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.