



Recombinant *Saccharomyces cerevisiae* Triosephosphate isomerase (TPI1)

Product Code	CSB-EP024102SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P00942
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	ARTFFVGGN FKLNGSKQSI KEIVERLNTA SIPENVEVVI CPPATYLDYS VSLVKKPQVT VQAQNAYLKA SGAFTGENSV DQIKDVGAKW VILGHSERRS YFHEDDKFIA DKTKFALGQG VGVILCIGET LEEKKAGKTL DVVERQLNAV LEEVKDWTNV VVAYEPVWAI GTGLAATPED AQDIHASIRK FLASKLGDKA ASELRILYGG SANGSNAVTF KDKADVDGFL VGGASLKPEF VDIINSRN
Source	E.coli
Target Names	TPI1
Protein Names	Recommended name: Triosephosphate isomerase Short name= TIM EC= 5.3.1.1 Alternative name(s): Triose-phosphate isomerase
Expression Region	2-248
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.