



Recombinant *Cricetulus griseus* Glutamine synthetase (GLUL)

Product Code	CSB-YP009553DXU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P04773
Product Type	Recombinant Protein
Immunogen Species	<i>Cricetulus griseus</i> (Chinese hamster) (<i>Cricetulus barabensis griseus</i>)
Purity	>85% (SDS-PAGE)
Sequence	<p>ATSASSHLN KGIKQMYMSL PQGEKVQAMY IWVDGTGEGL RCKTRTLDC</p> <p>PKCVEELPEW NFDGSSTFQS ESSNSDMYLS PVAMFRDPFR KEPNKLVFCE</p> <p>VFKYNQKPAE TNLRHTCKRI MDMVSNQHPW FGMEQEYTL</p> <p>GTDGHPFGWP SDGFPGPQGL YYCGVGADKA YRRDIMEAHY RACLYAGVKI</p> <p>TGTYAEVKHA QWEFQIGPCE GIRMGDHLWV ARFILHRVCK DFGVIATFDS</p> <p>KPIPGNWNGA GCHTNFSTKT MREENGLKHI KEAIEKLSKR HRYHIRAYDP</p> <p>KGGLDNARRL TGFHKTSNIN DFSAGVADRS ASIRIPRTVG QEKKGYFEAR</p> <p>CPSANCDPFA VTEAIVRTCL LNETGDQPFQ YKN</p>
Source	Yeast
Target Names	GLUL
Protein Names	<p>Recommended name: Glutamine synthetase Short name= GS EC= 6.3.1.2</p> <p>Alternative name(s): Glutamate decarboxylase EC= 4.1.1.15 Glutamate--ammonia ligase</p>
Expression Region	2-373
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>