



# Recombinant *Saccharomyces cerevisiae* Hexokinase-2 (HXK2)

<b>Product Code</b>	CSB-EP010469SVG
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P04807
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	VHLGPKKPQ ARKGSMADVP KELMQQIENF EKIFTVPTET LQAVTKHFIS ELEKGLSKKG GNIPMIPGWV MDFPTGKESG DFLAIDLGGT NLRVVLVKLG GDRTFDTTQS KYRLPDAMRT TQNPDELWEF IADSLKAFID EQFPQGISEP IPLGFTFSFP ASQNKINEGI LQRWTKGFDI PNIEHHDVVP MLQKQITKRN IPIEVVALIN DTTGTLVASY YTDPETKMGV IFGTGVNGAY YDVCSDIEKL QGKLSDDIPP SAPMAINCEY GSFENEHVVL PRTKYDITID EESPRPGQQT FEKMSSGYL GEILRLALMD MYKQGFIFKN QDLSKFDKPF VMDTSYPARI EEDPFENLED TDDLQNEFG INTTVQERKL IRRLSELIGA RAARLSVCGI AAICQKRGYK TGHIAADGSV YNRYPGFKEK AANALKDIYG WTQTSRDDYP IKIVPAEDGS GAGAAVIAAL AQKRIAEGKS VGIIGA
<b>Source</b>	E.coli
<b>Target Names</b>	HXK2
<b>Protein Names</b>	Recommended name: Hexokinase-2 EC= 2.7.1.1 Alternative name(s): Hexokinase PII Hexokinase-B
<b>Expression Region</b>	2-486
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.