



Recombinant Pig Chymotrypsin-like elastase family member 2A (CELA2A)

Product Code	CSB-EP005196PI-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P08419
Product Type	Recombinant Protein
Immunogen Species	Sus scrofa (Pig)
Purity	>85% (SDS-PAGE)
Sequence	VV GGEDARPSW PWQVSLQYDS SGQWRHTCGG TLVDQSWVLT AAHCSSSRT YRVVLGRHSL STNEPGSLAV KVKLVVHQD WNSNQLSNGN DIALKCLASP VSLTDKIQLG CLPAAGTILP NNYVCYVTGW GRLQTNGASP DILQQGQLLV VDYATCSKPG WVGSTVKTNM ICAGGDGIIS SCNGDSGGPL NCQGANGWQ VHGIVSFGSS LGCNYYHKPS VFTRVSNYID WINSVIANN
Source	E.coli
Target Names	CELA2A
Protein Names	Recommended name: Chymotrypsin-like elastase family member 2A EC=3.4.21.71 Alternative name(s): Elastase-2 Elastase-2A
Expression Region	29-269
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.