



Recombinant *Drosophila melanogaster* Nucleoside diphosphate kinase (awd)

Product Code	CSB-MP015889DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P08879
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	AANKERTFI MVKPDGVQRG LVGKIIERFE QKGFKLVALK FTWASKELLE KHYADLSARP FFPGLVNYMN SGPVVP MVWE GLNVVKTGRQ MLGATNPADS LPGTIRGDFC IQVGRNIIHG SDAVESAEKE IALWFNEKEL VTWTPAAKDW IYE
Source	Mammalian cell
Target Names	awd
Protein Names	Recommended name: Nucleoside diphosphate kinase Short name= NDK Short name= NDP kinase EC= 2.7.4.6 Alternative name(s): Abnormal wing disks protein Killer of prune protein
Expression Region	2-153
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.