



# Recombinant N-acetylneuraminase lyase (nanA)

<b>Product Code</b>	CSB-EP015994SZB-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P0A6L6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Shigella flexneri
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ATNLRGVMA ALLTPFDQQQ ALDKASLRRL VQFNIIQQGID GLYVGGSTGE AFVQSLSERE QVLEIVAEAA KGKIKLIAHV GCVSTAESQQ LAASAKRYGF DAVSAVTPFY YPFSFEEHCD HYRAIDSAD GLPMVVYNIP ALSGVKLTLD QINTLVTLPG VGALKQTS GD LYQMEQIRRE HPDLVLYNGY DEIFASGLLA GADGGIGSTY NIMGWRYQGI VKALKEGDIQ TAQKLQTECN KVIDLLIKTG VFRGLKTVLH YMDVVSVPLC RKPFGPVDEK YLPELKALAQ QLMQERG
<b>Source</b>	E.coli
<b>Target Names</b>	nanA
<b>Protein Names</b>	Recommended name: N-acetylneuraminase lyase EC= 4.1.3.3 Alternative name(s): N-acetylneuraminase pyruvate-lyase N-acetylneuraminic acid aldolase NALase Sialate lyase Sialic acid aldolase Sialic acid lyase
<b>Expression Region</b>	2-297
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.