



# Recombinant Ribosome-recycling factor (frr)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-BP014941SZB   |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | P0A808  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Shigella flexneri   |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | MISDIRKDAE VRMDKCVEAF KTQISKIRTG RASPSLLDGI VVEYYGTPTP<br>LRQLASVTVE DSRTLKINVF DRSMSPAVEK AIMASDLGLN PNSAGSDIRV<br>PLPPLTEERR KDLTKIVRGE AEQARVAVRN VRRDANDKVK ALLKDKEISE<br>DDRRSQDDV QKLTDAAIKK IEAALADKEA ELMQF   |
| <b>Source</b>            | Baculovirus   |
| <b>Target Names</b>      | frr   |
| <b>Protein Names</b>     | Recommended name: Ribosome-recycling factor Short name= RRF Alternative name(s): Ribosome-releasing factor  |
| <b>Expression Region</b> | 1-185   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |