



# Recombinant D-3-phosphoglycerate dehydrogenase (serA)

<b>Product Code</b>	CSB-BP017920EGX
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P0A9T1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AKVSLEKDK IKFLLVEGVH QKALESRAA GYTNIEFHKG ALDDEQLKES IRDAHFIGLR SRTHLTEDVI NAAEKLVAIG CFCIGTNQVD LDAAAKRGIP VFNAPFSNTR SVAELVIGEL LLLLRGVPEA NAKAHRGVWN KLAAGSFEAR GKKLGIIGYG HIGTQLGILA ESLGMYVYFY DIENKLPLGN ATQVQHLSDL LNMSDVVSLH VPENPSTKNM MGAKEISLMK PGSLINASR GTVVDIPALC DALASKHLAG AAIDVFPTEP ATNSDPFTSP LCEFDNVLLT PHIGGSTQEA QENIGLEVAG KLIKYSNNGS TLAVERNFEV SLPLHGGRRRL MHIHENRPGV LTALNKIFAE QGVNIAAQYL QTS AQMGYVV IDIEADEDVA EKALQAMKAI PGTIRARLLY
<b>Source</b>	Baculovirus
<b>Target Names</b>	serA
<b>Protein Names</b>	Recommended name: D-3-phosphoglycerate dehydrogenase Short name= PGDH EC= 1.1.1.95
<b>Expression Region</b>	2-410
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.