



# Recombinant Human Serpin B8 (SERPINB8)

<b>Product Code</b>	CSB-EP021076HU-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P50452
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDDLCEANGT FAISLFKILG EEDNSRVFF SPMSISSALA MVFMGAKGST AAQMSQALCL YKDGDHHRGF QSLLSEVNRT GTQYLLRTAN RLFGEKTCDF LPDFKEYCQK FYQAELEELS FAEDTEECRK HINDWVAEKT EGKISEVLDA GTVDPPLTKLV LVNAIYFKGK WNEQFDRKYT RGMLFKTNEE KKTVMQMMFKE AKFKMGYADE VHTQVLELPY VEEELSMVIL LPDDNTDLAV VEKALTYEKF KAWTNSEKLT KSKVQVFLPR LKLEESYDLE PFLRRLGMID AFDEAKADFS GMSTEKNVPL SKVAHKCFVE VNEEGTEAAA ATAVVRNSRC SRMEPRFCAD HPFLFFIRHH KTNCILFCGR FSSP
<b>Source</b>	E.coli
<b>Target Names</b>	SERPINB8
<b>Protein Names</b>	Recommended name: Serpin B8 Alternative name(s): Cytoplasmic antiproteinase 2 Short name= CAP-2 Short name= CAP2 Peptidase inhibitor 8 Short name= PI-8
<b>Expression Region</b>	1-374
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.