



# Recombinant Delta-aminolevulinic acid dehydratase (hemB)

<b>Product Code</b>	CSB-YP001558SZB
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P0ACB3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Shigella flexneri
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>TDLIQRPRR LRKSPALRAM FEETTLSLND LVLPIFVEEE IDDYKAVEAM PGVMRIPEKH LAREIERIAN AGIRSVMTFG ISHHTDETGS DAWREDGLVA RMSRICKQTV PEMIVMSDTC FCEYTSHGHC GVLCEHGVDN DATLENLGKQ AVVAAAAGAD FIAPSAAMDG QVQAIRQALD AAGFKDTAIM SYSTKFASSF YGPFREAAAGS ALKGDRKSYQ MNPMNRREAI RESLLDEAQQ ADCLMVKPKG AYLDIVRELR ERTELPIGAY QVSGEYAMIK FAALAGAIDE EKVVLESLGS IKRAGADLIF SYFALDLAEK KILR</p>
<b>Source</b>	Yeast
<b>Target Names</b>	hemB
<b>Protein Names</b>	Recommended name: Delta-aminolevulinic acid dehydratase Short name= ALAD Short name= ALADH EC= 4.2.1.24 Alternative name(s): Porphobilinogen synthase
<b>Expression Region</b>	2-324
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.