



Recombinant Human Palmitoyl-protein thioesterase 1 (PPT1)

Product Code	CSB-YP018587HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P50897
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	DPP APLPLVIWHG MGDSCCNPLS MGAIKKMVEK KIPGIYVLSL EIGKTLMEDV ENSFFLNVNS QVTTVCQALA KDPKLQQGYN AMGFSQGGQF LRAVAQRCPSPPMINLISVG GQHQGVFGLP RCPGESSHIC DFIRKTLNAG AYSKVVQERL VQAEYWHMPI KEDVYRNHSI FLADINQERG INESYKKNLM ALKKFVMVKF LNDSIVDPVD SEWFGFYRSG QAKETIPLQE TSLYTQDRLG LKEMDNAGQL VFLATEGDHL QLSEEWFYAH IIPFLG
Source	Yeast
Target Names	PPT1
Protein Names	Recommended name: Palmitoyl-protein thioesterase 1 Short name= PPT-1 EC= 3.1.2.22 Alternative name(s): Palmitoyl-protein hydrolase 1
Expression Region	28-306
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.