



Recombinant *Saccharomyces cerevisiae* 60S ribosomal protein L28 (RPL28)

Product Code	CSB-MP020217SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P02406
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	PSRFTKTRK HRGHVSAGKG RIGKHRKHPG GRGMAGGQHH HRINMDKYHP GYFGKVGMR Y FHKQQAHF WK PVLNLDKLWT LIPEDKRDQY LKSASKETAP VIDTLAAGYG KILGKGRIPN VPVIVKARFV SKLAE EKIRA AGGVVELIA
Source	Mammalian cell
Target Names	RPL28
Protein Names	Recommended name: 60S ribosomal protein L28 Alternative name(s): L27a L29 RP44 RP62 YL24
Expression Region	2-149
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.