



# Recombinant Papio anubis Protein Wnt-2 (WNT2)

<b>Product Code</b>	CSB-EP026133ERI
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A0M8S1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Papio anubis (Olive baboon)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SWWYM RATGGSSRVM CDNVPGLVSS QRQLCHRHPD VMRAISQGVA EWTAECQHGF RQHRWNCNTL DRDHSLFGRV LLRSSRESAF VYAISSAGVV FAITRACSQG EVKSCSCDPK KMGSAKDSKG IFDWGGCSDN IDYGIKFARA FVDAKERKGG DARALMNLHN NRAGRKAVKR FLKQECKCHG VSGSCTLRTC WLAMADFRKT GDYLWRKYNG AIQVVMNQDG TGFTVANERF KKPTKNDLVY FENSPDYCIR DREAGSLGTA GRVCNLTSRG MDSCEVMCCG RGYDTSHVTR MTKCGCKFHW CCAVRCQDCL EALDVHTCKA PKNADWTTPT
<b>Source</b>	E.coli
<b>Target Names</b>	WNT2
<b>Protein Names</b>	Recommended name: Protein Wnt-2
<b>Expression Region</b>	26-360
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.