



Recombinant Danio rerio RILP-like protein 1 (rilpl1)

Product Code	CSB-BP019713DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A0PJT0
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MEEQQLGAAL DKSAEELSVM DVYDIAAALG LELERVIERT GAELLSRLVP RVVRLVLELLE VLVSRSSSSP DTDELRLLELD RLRLRLRLERL EKEKKHKKEL ELVEDVWRGE AQDLLCQISQ LQEENKTLN NLSIRESPLT EEDIQKQEGM TERERQVMKK LKEVVDKQRD EIRAKDRELT LKNDDVEALQ QQMSRLMKIN QDVRHRVSVV EAQGKSLIQQ KVELEAAAQT QQQEVSSLRQ EVSRLKEKLLK EQSRSNEEEA QEPVGPSPSPA QEALCDEDLS TVDLKDPNRP RFTLQELRDV LHERNELKAK VFMLQEEIAY YRSEEQEEEN GPPLPDPSET LRTNPRSNFQ PESGIKRLFS FFSRDRSVSQ RRMMLNVEPV GDAVGSWTGK QEDVYTETAQ EALQHM
Source	Baculovirus
Target Names	rilpl1
Protein Names	Recommended name: RILP-like protein 1 Alternative name(s): Rab-interacting lysosomal-like protein 1
Expression Region	1-406
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.