



# Recombinant NAD-dependent malic enzyme (maeA)

<b>Product Code</b>	CSB-EP013633EJE
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A1AB75
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O1:K1 / APEC
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MEPKTKKQRS LYIPYAGPVL LEFPLLKGS AFSMEERRNF NLLGLLPEVV ETIEEQAERA WIQYQGFKTE IDKHIYLRNI QDNETLFYR LVNNHLDEMM PVIYTPTVGA ACERFSEIYR RSRGVFISYQ NRHNMDLILQ NVPNHNIKVI VVTGGERILG LGDQGIGGMG IPIGKLSLYT ACGGISPAYT LPVVLDVGTN NQQLLNPLY MGWRNPRITD DEYYEFVDEF IQAVKQRWPD VLLQFEDFAQ KNAMPLNRY RNEICSFNDD IQGTAAVTVG TLIAASRAAG GQLSEKKIVF LGAGSAGCGI AEMIIAQTR EGLSEEAARQ KVMVDRFGL LTKMKNLLP FQTKLVQKRE NLSDWTDSD VLSLLDVVRN VKPDILIGVS GQTGLFTEEI IREMHKHCPR PIVMPLSNPT SRVEATPQDI IAWTEGNALV ATGSPFNPVV WKDKIYPIAQ CNNAFIFPGI GLGVIASGAS RITDEMLMSA SETLAQYSPL VLNGEGLVLP ELKDIQKVSRAIAFAVGKMA QQQGVAVKTS AEALQQAIDD NFWHAEYRDY RRTSI
<b>Source</b>	E.coli
<b>Target Names</b>	maeA
<b>Protein Names</b>	Recommended name: NAD-dependent malic enzyme Short name= NAD-ME EC= 1.1.1.38
<b>Expression Region</b>	1-565
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.