



Recombinant *Pelobacter propionicus* Lipoyl synthase (lipA)

| | |
|--------------------------|---|
| Product Code | CSB-YP012927PYF |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | A1AMT5 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Pelobacter propionicus</i> (strain DSM 2379 / NBRC 103807 / OttBd1) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MDVQPTTMRG AVKTAHLGTD VSLATKPLPK PEWLRGKSAS TPDVERLVRI LRDNLRLHTVC EEASCPNLGE CFRKGTATFM IMGDVCTRHC PFCNVAHGSP HELAADEPVN LARAVELLKL SYVVITSVTR DDLPDGGAGH YGACVRALRD LKRSCLKVEIL TPDFRGAVAV AFEELRMNLP DVFNHNLETV PRLYPRVRPQ ADYHGSLDLL LRFREQFDHV PTKSGLMLGL GETEQEVRDV MEELRRHRCD MLTLGQYMRP SPHHLPEVERY VTPDEFERYR QFGLSIGFSH VESGPMVRSS YHADMQAREL MLVN |
| Source | Yeast |
| Target Names | lipA |
| Protein Names | Recommended name: Lipoyl synthase EC= 2.8.1.8 Alternative name(s): Lip-syn Short name= LS Lipoate synthase Lipoic acid synthase Sulfur insertion protein LipA |
| Expression Region | 1-314 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |