



Recombinant Human Nuclear receptor subfamily 2 group C member 2 (NR2C2)

Product Code	CSB-BP016052HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P49116
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>MTSPSPRIQI ISTDSAVASP QRIQIVTDQQ TGQKIQIVTA VDASGSPKQQ FILTSPDGAG TGKVLASPE TSSAKQLIFT TSDNLVPGRI QIVTDSASVE RLLGKTDVQR PQVVEYCVVC GDKASGRHYG AVSCEGCKGF FKRSVRKNLT YSCRSNQDCI INKHHRNRCQ FCRLKKCLEM GMKMESVQSE RKPFDVQREK PSNCAASTEK IYIRKDLRSP LIATPTFVAD KDGARQTGLL DPGMLVNIQQ PLIREDGTVL LATDSKAETS QGALGTLANV VTSLANLSES LNNGDTSEIQ PEDQSASEIT RAFDTLAKAL NTTDSSSSPS LADGIDTSGG GSIHVISRQDQ STPIIEVEGP LLSDTHVTFK LTMPSPMPEY LNVHYICESA SROLLFLSMHW ARSIPAFQAL GQDCNTSLVR ACWNEFLTLG LAQCAQVMSL STILAAIVNH LQNSIQEDKL SGDRIKQVME HIWKLQEFEN SMAKLDIDGY EYAYLKAVL FSPDHPGLTS TSQIEKFQEK AQMELQDYVQ KTYSEDYRL ARILVRLPAL RLMSSNITEE LFFTGLIGNV SIDSIIPIYL KMETAEYNGQ ITGASL</p>
Source	Baculovirus
Target Names	NR2C2
Protein Names	Recommended name: Nuclear receptor subfamily 2 group C member 2 Alternative name(s): Orphan nuclear receptor TAK1 Orphan nuclear receptor TR4 Testicular receptor 4
Expression Region	1-596
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.