



Recombinant Human 4-trimethylaminobutyraldehyde dehydrogenase (ALDH9A1)

Product Code	CSB-MP001582HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P49189
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	STGTFVVSQPLNYRGGARVEPADASGTEKAFEPATGRVIATFTCSGEKEVNLA VQNAKA AFKIWSQKSGMERCILLEAARIIREREDEIATMECINNGKSIFEARLDI DISWQCLEY YAGLAAS MAGEHIQLP GGSFGYTRREPLGVCV GIGAWNYPFQIA SWKSAPALACGNAMVFKPSPFTPVSALLLAEIYSEAGVPPGLFN VVQGGAA TG QFLCQHPDVAKVSFTG SVPTGMKIMEMSAKG IKPVTLELGGKSPLIIFSDCDMN NAVKGALMANFLTQ GQVCCNGTRV FVQKEILDKFTEE VVKQTQRIKIGDPLLE DTRMGPLINRPHLERVLGFVKVAKEQGAKVLCGGDIYVPEDPKLKDGY YMRP CVLTNCRDDMT CVKEEIFGPVMSILSFDTEAEVLERANDTT FGLAAGVFTRDIQ RAHRVVAELQAGTCFINN YNVSPVELPFGG YKKSFGRENGRVTIEY SQLKT VCVEMGDVESAF
Source	Mammalian cell
Target Names	ALDH9A1
Protein Names	Recommended name: 4-trimethylaminobutyraldehyde dehydrogenase Short name= TMABADH EC= 1.2.1.47 Alternative name(s): Aldehyde dehydrogenase E3 isozyme Aldehyde dehydrogenase family 9 member A1 EC= 1.2.1.3 Gamma-aminobut
Expression Region	2-494
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein belongs to the aldehyde dehydrogenase family of proteins. It has a high activity for oxidation of gamma-aminobutyraldehyde and other amino aldehydes. The enzyme catalyzes the dehydrogenation of gamma-aminobutyraldehyde to gamma-aminobutyric acid (GABA). This isozyme is a tetramer of identical 54-kD subunits.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.