



Recombinant Danio rerio Ataxin-7-like protein 3 (atxn7l3)

Product Code	CSB-EP002448DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A1L209
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MKMEDMSLSG LDNTKLEALA HDVYSDLVED ACLGLCFEVH RAVKQGYFFL DETDQESMKD FEIVDQPGVD IFGQVYNQWK NKECVCPCNCS RSIAASRFAP HLEKCLGMGR NSSRIANRRI ASSNNTSKSE SDQEDNDDIN DNDWSYGSEK KAKKRKSEKN PNSPRRSKSL KHKNGELSGS VNPDMYKYN Y SSGISYETLG PEELRSILTT QCGVVSEHTK KMCTRSQRCP QHTDEQRRAV RVFLLGPSAS TLPDADTMLE NEAYEPPDGQ LIMSRLHWDA SSDISPSDSA SSKASTNNS SKRPKKKKPS TSLTPAGER DKAQERDRIA GSGSSGSSSQ NALGLSSRKK RPKLAVPPAP SIYDDLN
Source	E.coli
Target Names	atxn7l3
Protein Names	Recommended name: Ataxin-7-like protein 3 Alternative name(s): SAGA-associated factor 11 homolog
Expression Region	1-367
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.