



Recombinant *Shewanella amazonensis* Phosphatidylserine decarboxylase proenzyme (psd)

Product Code	CSB-MP018033STK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A1SA30
Product Type	Recombinant Protein
Immunogen Species	<i>Shewanella amazonensis</i> (strain ATCC BAA-1098 / SB2B)
Purity	>85% (SDS-PAGE)
Sequence	MDKVKIALQY LMPKHLVSRI VGKFAAAEAG FVTTAFIKWF IKQYGINMSE ALHSNPEAYK TFNDFFTREL KPGLRPIDQA EDIMVHPVDG AVSQLGPIED GRIFQAKGHH YSALALLGGQ ADDAARFEEG DFATIYLAPK DYHRIHMPIK GTLSKMTYVP GELFSVNPLT ARNVPGLFAR NERVVAIFET DKGPLAMVLV GATIVASJET VWAGTVTPPT GKKVFTWDYP TEGPEALTLE KGAEMGRFKL G
Source	Mammalian cell
Target Names	psd
Protein Names	Recommended name: Phosphatidylserine decarboxylase proenzyme EC=4.1.1.65 Cleaved into the following 2 chains: 1. Phosphatidylserine decarboxylase alpha chain 2. Phosphatidylserine decarboxylase beta chain
Expression Region	1-251
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.