



# Recombinant Human LIM/homeobox protein Lhx2 (LHX2)

<b>Product Code</b>	CSB-EP012920HU-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P50458
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLFHSLSGPE VHGVIDEMDR RAKSEAP AIS SAIDRGDTET TMPSISSDRA ALCAGCGGKI SDRYYLLAVD KQWHMRCLKC CECKLNLESE LTCFSKDGSI YCKEDYYRRF SVQRCARCHL GISASEMVMR ARDLVYHLNC FTCTTCNKML TTGDHFGMKD SLVYCRHLFE ALLQGEYPAH FNHADVAAAA AAAAAAKSAG LGAAGANPLG LPYYNGVGT V QKGRPRKRKS PGPGADLAAY NAALSCNEND AEHLDRDQPY PSSQKTKRMR TSFKHHQLRT MKSYFAINHN PDAKDLKQLA QKTGLTKRVL QVWFQNARAK FRRNLLRQEN TGVDKSTDAALQTGTPSGPA SELSNASLSP SSTPTTLTDL TSPTLPTVTS VLTSVPGNLE GHEPHSPSQT TLTNLF
<b>Source</b>	E.coli
<b>Target Names</b>	LHX2
<b>Protein Names</b>	Recommended name: LIM/homeobox protein Lhx2 Short name= Homeobox protein LH-2 Short name= LIM homeobox protein 2
<b>Expression Region</b>	1-406
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	This gene encodes a protein belonging to a large protein family, members of which carry the LIM domain, a unique cysteine-rich zinc-binding domain. The encoded protein may function as a transcriptional regulator. The protein can recapitulate or rescue phenotypes in Drosophila caused by a related protein, suggesting conservation of function during evolution.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.