



# Recombinant Human Granzyme M (GZMM)

<b>Product Code</b>	CSB-MP010085HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P51124
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	IIGGR EVIPHSRPYM ASLQRNGSHL CGGVLVHPKW VLTAAHCLAQ RMAQLRLVLG LHTLDSPGLT FHIKAAIQHP RYKVPVALEN DLALLQLDGK VKPSRTIRPL ALPSKRQVVA AGTRCSMAGW GLTHQGGRLS RVLRELDLQV LDTRMCNNSR FWNGSLSPSM VCLAADSKDQ APCKGDSGGP LVCGKGRVLA RVLSFSSRVC TDFIKPPVAT AVAPYVSWIR KVTGRSA
<b>Source</b>	Mammalian cell
<b>Target Names</b>	GZMM
<b>Protein Names</b>	Recommended name: Granzyme M EC= 3.4.21.- Alternative name(s): Met-1 serine protease Short name= Hu-Met-1 Met-ase Natural killer cell granular protease
<b>Expression Region</b>	26-257
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	Human natural killer (NK) cells and activated lymphocytes express and store a distinct subset of neutral serine proteases together with proteoglycans and other immune effector molecules in large cytoplasmic granules. These serine proteases are collectively termed granzymes and include 4 distinct gene products: granzyme A, granzyme B, granzyme H, and Met-ase, also known as granzyme M.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.