



Recombinant Echinops telfairi Protein Wnt-2 (WNT2)

Product Code	CSB-YP026133EJK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A1X153
Product Type	Recombinant Protein
Immunogen Species	Echinops telfairi (Lesser hedgehog tenrec)
Purity	>85% (SDS-PAGE)
Sequence	SWWYM RATGTSRVMC DNVPLVSRQ RQLCHRHPDV MRAIGLVAE WTAECQHQRFR QHRWNCDTLD RDHSLFGRLL LRSSRESAFV YAISSAGVVF AITRACSQGE LKSCSCDPKK KGTAKDSRGT FDWGGCSDNI DYGVKFARAF VDAKEKKGKD ARALMNLHNN RAGRKAVKRF LKQECKCHGV SGSCTLRTCW LAMADFRKTG DYLRWRYNGA IQVVMNQDGT GFTVANKRFK KPTKNLTVYF ENSPDYCIRD RDAGSPGTAG RVCNLTSRGM DSCEVMCCGR GYDTSRVTRM TKCECKFHWC CAVRCQDCLE ALDVHTCKAP KSADWAVPT
Source	Yeast
Target Names	WNT2
Protein Names	Recommended name: Protein Wnt-2
Expression Region	26-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.