



Recombinant Human Zinc finger protein 75D (ZNF75D)

Product Code	CSB-MP027016HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P51815
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAMRELNADS CSSPQMGAMW ETSGSVKENS SQSKKYSTKI ENLGPESACR HFWSFRYHEA TGPLETISQL QKLCHQWLRP EIHSKEQILE MLVLEQFLSI LPKETQNWVQ KHHPQNVKQA LVLVEFLQRE PDGTKNEVTA HELKKEAVLL GGTAVAPGFK WKPAEPQPMG VFQKEYWNTY RVLQEQLGWN THKETQPVYE RAVHDQQMLA LSEQKRIKHW KMASKLILPE SLSLLTFEDV AVYFSEEEWQ LLNPLEKTLY NDVMQDIYET VISLGLKLN DTGNDHPISV STSEIQTSGC EVSKKTRMKI AQKTMGRENP GDTHSVQKWH RAFPRKKRKK PATCKQELPK LMDLHGKGPT GEKPFKCQEC GKSFRVSSDL IKHHRIHTGE KPYKCQQCDR RFRWSSDLNK HFMTHQGIKP YRCSWCGKSF SHNTNLHTHQ RIHTGEKPFK CDECGRFIQ NSHLIKHQRT HTGEQPYTCS LCKRNFSRRS SLLRHQKLHR RREACLVSPN
Source	Mammalian cell
Target Names	ZNF75D
Protein Names	Recommended name: Zinc finger protein 75D Alternative name(s): Zinc finger protein 75 Zinc finger protein 82
Expression Region	1-510
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.