



Recombinant Human Heterogeneous nuclear ribonucleoprotein A3 (HNRNPA3)

Product Code	CSB-MP010603HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P51991
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEVKPPPGRP QPDSGRRRRR RGEEGHDPKE PEQLRKLFIG GLSFETTTDDSLREHFEKWGT LTDCVVMRDP QTKRSRGGFGF VTYSCEEVD AAMCARPHKV DGRVVEPKRA VSREDSVKPG AHLTVKKIFV GGIKEDTEEYNLRDYFEKYG KIETIEVMED RQSGKKRGFA FVTFDDHDTV DKIVVQKYHT INGHNCEVKK ALSKQEMQSA GSQRGRGGGS GNFMGRGGNF GGGGGNFGRG GNFGGRGGYG GGGGSRGSY GGDGGYNGF GGDGGNYGGG PGYSSRGGYG GGGPGYGNQG GYGGGGGYD GYNEGGNFGG GNYGGGGNYN DFGNYSGQQQ SNYGPMKGGG FGGRSSGSPY GGGYGSGGGG GYGSRRF
Source	Mammalian cell
Target Names	HNRNPA3
Protein Names	Recommended name: Heterogeneous nuclear ribonucleoprotein A3 Short name= hnRNP A3
Expression Region	1-378
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.