



Recombinant Danio rerio HORMA domain-containing protein 1 (hormad1)

Product Code	CSB-EP010644DIL-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A2BF66
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MTCEQKSRPV QSCQMLPDVV SNEQQSLVLV KKLLAIAISS ITYLRGLFSE KAYGRKYVGE LKVYILREDS CPGAQQIVHW LQGCFDALQR RYLRMVLLSI YCDPDNPQKV TECYQFKIKY TEKGPQMDFE SKNGQSLTKM ACDNTQKSSM LLVRKLYMLM QNLGPLPDDV CLNMKLLYYD EVTPQEYQPP GFKEDDNGTL MFEREPVNL MGEVVTPFHS IKMNVTEKK RIEPFEDDVE VCVSTKWSLK ISEDGMMSET SVQQECMTKE NVITDAGIEY SETQETQEQP HRHTKEDFST NPKMDNLVKK TADLKVDARK TRSGRIFEPQ ISQLEFPLSQ DPQPSVPKRR KVSVPK
Source	E.coli
Target Names	hormad1
Protein Names	Recommended name: HORMA domain-containing protein 1
Expression Region	1-356
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.