



Recombinant Danio rerio Coiled-coil domain-containing protein 125 (ccdc125)

Product Code	CSB-MP004645DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A2BGP7
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MGARSRRNHG HLPQMKTRGR FLDFLSHSRP IKEGGERSNG PFSWTSCRAM YSVFREELEA ERVARWEKRQ SESSTEDTSE ELQRRRLQEV EEVELLRTTEL EVTHRHLEGG HEALRILQGQ AILDKATCHT KMLLQKSEER TKALEKEVNA LQWEITFNQV QFKNVENSWS LKYERVLAEN EALKKGLEEK MKEHQKQRTE NASLSQKCLE LLSMLSAKER RDFQRTQPSC SLRTDGSAL LAVYGACQCN SNGGEPSCA RSAAASRKQV LQLKQELEQQ QKRKDEAYVM MDAFRIAFAQ QLRRVGENVL RQAETDRHQT HNQRHEKGGKQ WSLTVGERLK KILPTTSDGR IPSDSSETLH MLLDLLNDKE EALAHQRKVS YMLARNTENL EKRLLMQLEE LDVSHEESKT KTVASEDYEW SRDCRCCADS CGSQHLSVPD DKQTYKANQP PSENRPETKP DNKDA
Source	Mammalian cell
Target Names	ccdc125
Protein Names	Recommended name: Coiled-coil domain-containing protein 125
Expression Region	1-465
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.