



Recombinant Danio rerio Cytosolic Fe-S cluster assembly factor narfl (narfl)

Product Code	CSB-MP015455DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A2RRV9
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MAVIFRSISG CITSAESVLI TQQSHEELYR VLRHMKQVSS TEQKVVVSV SPQSRASLAA HYGIGSSEVA RKLTSFLKHL GVHHVFDTAF SRSFSLIESQ REFLQRFSQR EADKKALPML ASACPGWICY AEKTHGEFIL PYISTTRSPQ QIMGSLVKGY FASQKGVSPQ MIYHVTVMPC YDKKLEASRP DFYLSEHETR EVDCVITSGE VLKMLEEEKV SLRDVQPAPL DTMFSNVCGE ELLGHAGSGS GGYLHHIYKH AAKQLFGVDV DELTYKTMKN KDFQEVTLEK DGQVLLKFAA VYGFNRNIQNL VQKLKRGKSP YHFVEVMACP SGCLNGGGQL KPSADQSNKE LLQQVEEVYR SEHPSVPEDD SQVAELYQSW LESVGEEKAR QLLHTQYHAV EKTANGLSIK W
Source	Mammalian cell
Target Names	narfl
Protein Names	Recommended name: Cytosolic Fe-S cluster assembly factor narfl Alternative name(s): Nuclear prelamin A recognition factor-like protein
Expression Region	1-411
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.