



Recombinant *Oryza sativa* subsp. *indica* Autophagy protein 5 (ATG5)

Product Code	CSB-EP002293OFF-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A2X052
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MAAQRDDEAG WSAEAARRVW GGAVPLQVHL HDADVTTLPP PPPFLTGLGPR IGYLPLLVI IKAHFSSTLP PGIDTVWFEY KGLPLKWIYIP IGVLYDLLCA DPERPWNLTV HFRGYPSEIL TPCDGEDSVK WSYMNSLKEA AFIITGNSKN VMNMSQADQG ALWQSVMKGN LDGYMNISTR LKLGPFEEDC LVRTSSVEGQ QGSDEPESPG SGKPCRVVPR LYVRSVQEDL YDLEDALPVG DWESISYINR PFEVRREEGR SYITLEHALK TLLPEFFSSK ASRIPDDSET APQAPDSAPN DDSDVTPRSC EKLESSASSS PQEANVANKG KIVKLVRVQG IEVDMDIPFL WVANNLNKPE CYLHICVYVG TRKREPKDGR
Source	E.coli
Target Names	ATG5
Protein Names	Recommended name: Autophagy protein 5
Expression Region	1-380
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.