



Recombinant Danio rerio Oxidoreductase NAD-binding domain-containing protein 1 (oxnad1)

Product Code	CSB-YP017311DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A3KP77
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MTSRRRTDHL ERTASVHRQM ELFSARVCDI ISESDTVKRL RLEVAHPDFS FRAGQWVDF IPGVDTVGGF SICSSPGLLK REGAIELAVK YARHPPAHWI HTECSVDSQV AVRVGGNFYF DPQPSNPVVD LLLVAGGVGI NPLYSIILLHA ADLHRHSH RYTPGHSHLC YSAKNTTELL FKDTIIDICH ERPDKFSCHF HVTQQSSDIE PQLQPYTIRG RISAEELQRY VDPERTLCYL CGPPPMIEKV SSDLQSTGLP EDRILFEKWW
Source	Yeast
Target Names	oxnad1
Protein Names	Recommended name: Oxidoreductase NAD-binding domain-containing protein 1 EC= 1.-.-.
Expression Region	1-270
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.