



Recombinant Mouse Transcription factor Maf (Maf)

Product Code	CSB-YP013312MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P54843
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MASELAMNNS DLPTSPLAME YVNDFDLMKF EVKKEPVETD RIISQCGRLI AGGSLSTPM STPCSSVPPS PSFSAPSPGS GSEQKAHLED YYWMTGYPQQ LNPEALGFSP EDAVEALISN SHQLQGGFDG YARGAQQQLAA AAGAGAGASL GGSGEEMGPA AAVVSAVIAA AAAQSGAAPH YHHHHHHAAG HHHHPTAGAP GAAGGASASA SGAGGAGGGG PASAGGGGGG GGGGGTAGAG GALHPHHAAG GLHFDDRFS EQLVTMSVRE LNRQLRGVSK EEVIRLKQKR RTLKNRGYAQ SCRFRVQQR HVLESEKNQL LQQVDHLKQE ISRLVRERDA YKEYEKLVS NGFRENGSSS DNPSSPEFFM
Source	Yeast
Target Names	Maf
Protein Names	Recommended name: Transcription factor Maf Alternative name(s): Proto-oncogene c-Maf V-maf musculoaponeurotic fibrosarcoma oncogene homolog
Expression Region	1-370
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.