



# Recombinant *Oryza sativa* subsp. *indica* Alcohol dehydrogenase class-3 (ADHIII)

<b>Product Code</b>	CSB-MP001357OFF
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A2XAZ3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MASSTQGQVI TCKAAVAWEA NKPMTIEDVQ VAPPQAGEVR VKILFTALCH TDHYTWSGKD PEGLFPCILG HEAAGIVESV GEGVTEVQPG DHVIPCYQAE CRECKFCKSG KTNLCGKVRA ATGVGVMMND RKSRSFINGK PIYHFMGTST FSQYTVVHDV SVAKINPQAP LDKVCLLGG VSTGLGAVWN TAKVEAGSIV AIFGLGTVGL AVAEGA KASAG ASRIIGIDID SKKFDVAKNF GVTEFVNP KD HDKPIQQVIV DLT DGGVDYS FECIGNVSVM RSALECCHKG WGTSVIVGVA ASGQEISTRP FQLVTGRVWK GTA FGGFKSR SQVPWLVEKY LNKEIKVDEY VTHSMNLTDI NKA FDLLHEG GCLRCVLATD K
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ADHIII
<b>Protein Names</b>	Recommended name: Alcohol dehydrogenase class-3 EC= 1.1.1.1 Alternative name(s): Alcohol dehydrogenase class-III Glutathione-dependent formaldehyde dehydrogenase Short name= FALDH Short name= FDH Short name= GSH-FDH
<b>Expression Region</b>	1-381
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.