



Recombinant Danio rerio 2-oxoglutarate and iron-dependent oxygenase domain-containing protein 2 (ogfod2)

Product Code	CSB-MP016310DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A3KGZ2
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	≥85% (SDS-PAGE)
Sequence	MERFYTCSCF FTDNIFLECY KLHVRVSEN QFRKDYQNIL RSLGCESESQ FRDVIGKIQA EIERRQNHKL KSTERAARIK EIYTPLHQHV YHLQESFLAP ELLEMVKYCA SSEANVQGLL KLIQTEAASR VFRFQVFRKE FCKDLLEELE HFEQSDAPKG RPNTMNNYGI VLNELGFDEG FITPLREVYL RPLTALLYSD CGGNCLDSHK AFVVKYDMHE DLNLSYHYDN SEVTNLNVLG KDFTEGNLFF GDMRQVPLSE TECVEVEHRV TEGLLHRGQH MHGALSISSG TRWNLIWMR ASRQRNKLCP MCGKRPTLVE SDGFSDGFTM DSDGDARANV SCSLT
Source	Mammalian cell
Target Names	ogfod2
Protein Names	Recommended name: 2-oxoglutarate and iron-dependent oxygenase domain-containing protein 2 EC= 1.14.11.-
Expression Region	1-345
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.