



Recombinant Macaca mulatta (Rhesus macaque) Prolactin (PRL)

Product Code	CSB-MP018724MOW
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P55151
Product Type	Recombinant Protein
Immunogen Species	Macaca mulatta (Rhesus macaque)
Purity	>85% (SDS-PAGE)
Sequence	LP ICPGGAARCQ VTLRDLFDRA VVLSHYIHNL SSEMFSEFDK RYTHGRGFIT RAINSCHTSS LPTPEDKEQA QQMNQKDFLS LIVSILRSWN EPLYHLVTEV RGMEEAPEAI LSKAVEIEEQ TKRLLGEMEL IVSQVHPETK ENEIYPVWTG LPSLQMADEE SRLSAYYNLL HCLRRDSHKI DNYLKLLKCR IHHNNC
Source	Mammalian cell
Target Names	PRL
Protein Names	Recommended name: Prolactin Short name= PRL
Expression Region	29-227
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.