



ABI2 Antibody

Product Code	CSB-PA001102GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYB9
Immunogen	Human ABI2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	abl interactor 2;ABI2;ABI-2;ABI2B;AIP-1;AbIBP3;SSH3BP2;argBP1A;argBP1B ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ABI2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.