



ACTB Antibody

Product Code	CSB-PA001207GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60709
Immunogen	Human ACTB
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat,Zebrafish
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	actin, beta;ACTB;PS1TP5BP1 ;beta-actin; ??-actin
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ACTB
Target Details	This gene encodes one of six different actin proteins. Actins are highly conserved proteins that are involved in cell motility, structure, and integrity. This actin is a major constituent of the contractile apparatus and one of the two nonmuscle cytoskeletal actins.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.