



# Recombinant Human Caspase-6 (CASP6), partial

<b>Product Code</b>	CSB-MP004551HU
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P55212
<b>Storage Buffer</b>	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AFYKREM FDPAEKYKMD HRRRGIALIF NHERFFWHLT LPERRGTCAD RDNLTRRFSD LGFEVKCFND LKAEELLLKI HEVSTVSHAD ADCFVCVFLS HGEGNHIYAY DAKIEIQTLT GLFKGDKCHS LVGKPKIFII QACRGNQHDV PVIPLDVVD
<b>Source</b>	Mammalian cell
<b>Target Names</b>	CASP6
<b>Protein Names</b>	Recommended name: Caspase-6 Short name= CASP-6 EC= 3.4.22.59 Alternative name(s): Apoptotic protease Mch-2 Cleaved into the following 2 chains: 1. Caspase-6 subunit p18 2. Caspase-6 subunit p11
<b>Expression Region</b>	24-179
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.