



ACTR3B Antibody

Product Code	CSB-PA001250GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P1U1
Immunogen	Human ACTR3B
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	ARP3 actin-related protein 3 homolog B (yeast);ACTR3B;ARP11;ARP3BETA;DKFZp686O24114 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ACTR3B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.