



# ADAT2 Antibody

<b>Product Code</b>	CSB-PA001328GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q7Z6V5
<b>Immunogen</b>	Human ADAT2
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB,IHC
<b>Storage Buffer</b>	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity purified
<b>Isotype</b>	IgG
<b>Alias</b>	adenosine deaminase, tRNA-specific 2, TAD2 homolog (S. cerevisiae);ADAT2;DEADC1;DKFZp686L1118;FLJ44213;TAD2;dJ20N2;dJ20N2.1 ;
<b>Product Type</b>	Purified Rabbit Anti human PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	ADAT2
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.