



# ADCY1 Antibody

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA001336GA01HU  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| <b>Uniprot No.</b>         | Q08828  |
| <b>Immunogen</b>           | Human ADCY1   |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Human,Mouse,Rat   |
| <b>Tested Applications</b> | ELISA,WB  |
| <b>Storage Buffer</b>      | 0.1M NaHCO <sub>3</sub> ,0.1M Glycine, 0.02% Sodium Azide, 50% Glycerol, pH 7.3.<br>-20°C, Avoid freeze / thaw cycles.                      |
| <b>Purification Method</b> | Antigen Affinity purified   |
| <b>Isotype</b>             | IgG   |
| <b>Alias</b>               | adenylate cyclase 1 (brain);ADCY1;AC1 ;   |
| <b>Product Type</b>        | Purified Rabbit Anti Human,Mouse PolyClonal Antibody  |
| <b>Immunogen Species</b>   | Homo sapiens (Human)  |
| <b>Target Names</b>        | ADCY1   |
| <b>Target Details</b>      | This gene encodes a form of adenylate cyclase expressed in brain. A similar protein in mouse is involved in pattern formation of the brain. |